

ATCC Medium 1395: Middlebrook 7H9 Broth w/ ADC Enrichment

Base Medium

Middlebrook 7H9 Broth.....	4.7 g
Glycerol.....	2.0 mL
DI Water.....	900.0 mL
ADC Enrichment (see below)	100.0 mL

Add 4.7 grams of dehydrated Middlebrook 7H9 medium to 900 mL of water containing 2.0 mL of glycerol.

Mix well and autoclave at 121°C. Cool to 50-55°C.

Aseptically add filter sterilized ADC Enrichment.

ADC Enrichment

Bovine Serum Albumin, Frac V.....	5.0 g
Dextrose.....	2.0 g
Catalase(beef).....	3.0 mg
DI Water.....	100.0 mL

Filter sterilize.