

**ATCC Medium: 1452 Mueller Hinton Medium with 10% Rabbit Serum**

**Broth**

Mueller Hinton Broth (BD 211443).....22.0 g  
DI Water.....900 ml

Autoclave at 121°C and let cool. Aseptically add 100 ml of sterile rabbit serum. Dispense into appropriate vessel.

**For Solid Media**

Mueller Hinton Agar (BD 225250).....38.0 g  
DI Water.....900 ml

Autoclave at 121°C and let cool to 48°C. Aseptically add 100 ml of sterile rabbit serum. Dispense into appropriate vessel.