

ATCC medium: 587 YPSC medium

Yeast extract.....1.0 g
Peptone.....1.0 g
Sodium acetate . 3H₂O.....0.5 g
L-Cysteine . HCl.....0.05 g
MgSO₄ . 7H₂O0.25 g
CaCl₂ . 2H₂O0.25 g
Distilled water.....1.0 L

Adjust pH to 7.5 with 10 M NaOH. For semisolid and solid media add 6.0 g or 15.0 g agar respectively per liter of the medium.
Autoclave at 121C for 15 minutes.